

Children's beds



Eleganza Smart Junior, Tom, Mimi

Children's beds

ELEGANZA SMART JUNIOR

- ACTIVE PATIENT MOBILISATION SUPPORT
- PERFECT PROTECTION AGAINST PATIENT FALLS
- SAFE TRANSFER OF PATIENT



The siderails have smaller gaps between the rails in comparison with standard beds.

FOR PATIENT AND NURSE

- The bed is fully electrically adjustable. The actuators can adjust the position of the 4-section mattress platform made of metal segments and of the bed height.
- Auto-regression of the backrest enables the enlargement of the pelvic area by 11 cm.
- The backup battery with an intelligent charger ensures the bed's safe operation.
- Unique GO safety concept prevents accidental bed positioning.

TECHNICAL PARAMETERS

External dimensions	92,5×183,5 cm
Mattress platform	80×170 cm
Mattress platform height adjustment	38,5 – 76,5 cm
Backrest tilt	70°
Thighrest tilt	45°
Backrest auto-regression	11 cm
TR/ATR tilt	±16°
Safe working load	180 kg
Weight (depending on the equipment)	120 kg
Clearance height	14,5 cm

MIMI



- The Mimi baby bed with height adjustment and bed surface tilt is suitable for the hospital neonatal wards
- Height adjustment and tilting to the Trendelenburg and ATR position ($\pm 14^\circ$)
- The bed consists of a transparent plastic removable tub
- Rooming-In is mobile: 4×65 mm diameter double wheels, 2 have separate brakes
- Dimensions: 58 cm width, 96 cm length, height with handrail 101–131 cm
- A washable mattress is supplied as an optional accessory

TOM



- Children's bed in variants of Tom 140 and Tom 115 with mechanical height adjustment of the bed surface and TR/ATR tilting position for standard children's wards
- Tom 140 allows positioning by a mechanical screw mechanism, dimensions of bed surface 70×136 cm (height without mattress 78–109 cm), outer dimensions 83×165 cm
- Tom 115 without back section positioning, bed surface dimensions 55×115 cm (height without mattress 78–109 cm), outer dimensions 69×144 cm
- The bed has an integrated lowering side rail adjustable to 3 positions. Adjustability is secured against accidental dropping or child manipulation.
- The Tom bed is mobile: 4 double 100 mm diameter wheels with separate brakes, and protective wheels in the bed corners and at head end
- 2 infusion stand holders, diagonally placed in the bed corners
- A washable mattress is supplied as an optional accessory